



Creating balance in healthcare



## Timo & Cleo Shower-toilet chair



# Timo Shower-toilet chair



5400 6020  
**Timo Shower chair**

- 2 independently foldable armrests.
- Upholstered seat.
- Snaplock removable seat.
- Seat 46 x 45 x 50 cm /  
18.1 x 17.7 x 19.7 inch (w x d x h).
- Maximum load 135 kg / 300 lbs.

5400 6021  
**Timo Shower-toilet chair**

- 2 independently foldable armrests.
- Upholstered seat with toilet opening.
- Universal rail for bedpan and / or toilet bucket.
- Snaplock removable seat.
- Seat 46 x 45 x 50 cm /  
18.1 x 17.7 x 19.7 inch (w x d x h).
- Maximum load 135 kg / 300 lbs.



5400 6025  
**Timo Shower chair mobile**

- 2 independently foldable armrests.
- Upholstered seat.
- 4 double braked castors.
- Snaplock removable seat.
- Seat 48 x 49 x 52 cm /  
18.9 x 19.3 x 20.5 inch (w x d x h).
- Maximum load 135 kg / 300 lbs.





Suitable for people with a mild to moderate disability



5400 6030

### **Timo Shower-toilet chair mobile**

- 2 independently foldable armrests.
- 4 double braked castors.
- Upholstered seat with toilet opening.
- Universal rail for bedpan and / or toilet bucket.
- Snaplock removable seat.
- Seat 48 x 49 x 54 cm /  
18.9 x 19.3 x 21.3 inch (w x d x h).
- Maximum load 135 kg / 300 lbs.



5400 6035

### **Timo Shower-toilet chair mobile + footrests**

- 2 independently foldable armrests.
- 4 double braked castors.
- Upholstered seat with toilet opening.
- Universal rail for bedpan and / or toilet bucket.
- Footrests foldable, swing-out and removable.
- Snaplock removable seat.
- Seat 48 x 49 x 54 cm /  
18.9 x 19.3 x 21.3 inch (w x d x h).
- Maximum load 135 kg / 300 lbs.

# Cleo Shower-toilet chair



Suitable for people with a mild to moderate disability



5400 6040

## Cleo Shower chair self-propelled

- 2 independently foldable armrests.
- Upholstered seat.
- Footrests foldable, swing-out and removable.
- 2 plastic 24 inch wheels with massive tires.
- 2 non-braked castors.
- Snaplock removable seat.
- Seat 46 x 45 x 52 cm /  
18.1 x 17.7 x 20.5 inch (w x d x h).
- Maximum load 135 kg / 300 lbs.



5400 6050

## Cleo Shower-toilet chair self-propelled

- 2 independently foldable armrests.
- Upholstered seat with toilet opening.
- Universal rail for bedpan and / or toilet bucket.
- Footrests foldable, swing-out and removable.
- 2 plastic 24 inch wheels with massive tires.
- 2 non-braked castors.
- Snaplock removable seat.
- Seat 46 x 45 x 54 cm /  
18.1 x 17.7 x 21.3 inch (w x d x h).
- Maximum load 135 kg / 300 lbs.

For more information, visit our website: [www.lopital.com](http://www.lopital.com)

Lopital Nederland B.V., Laarakkerweg 9, 5061 JR Oisterwijk, Postbus 56, 5060 AB Oisterwijk  
Tel +31 (0)13 5239300, Fax +31 (0)13 5239301, E-mail [info@lopital.nl](mailto:info@lopital.nl), [www.lopital.nl](http://www.lopital.nl)

